

MammoStar™ Tissue Markers

Ultrasound-Guided Wire Localization

Clinical History

51 year-old female presented for diagnostic breast imaging to evaluate a 1 week palpable mass in the right breast.

Findings

Mammogram and breast ultrasound revealed a highly suspicious 1.9cm mass in the 10:00 7 CFN-B position of the right breast. Ultrasound-Guided core biopsy was recommended.

Procedure & Pathology

Ultrasound guided vacuum-assisted Mammotome® EX biopsy was performed. A MammoStar™ was deployed at the biopsy site adjacent to the residual mass (*figure 1*). Pathology results revealed an infiltrating ductal carcinoma, grade III.

MRI Imaging


A preoperative staging breast MRI was performed. The enhancing residual malignant mass is identified adjacent to an oval signal void, which corresponds to the MammoStar™ at 8 days after deployment (*figure 2*).

Ultrasound-Guided Wire Localization


Ultrasound-guided wire localization for definitive surgery and sentinel lymph node evaluation was performed. Images of MammoStar™ at 12 days (*figure 3*) reveal hyperechoic convex lines (beta-glucan) adjacent to the residual malignant mass. The post-procedure mammogram (*figure 4*) confirmed accuracy of the wire placement.

Courtesy


Kimberly C. Hutcherson, M.D.
 North Metropolitan Radiology Associates, LLP
 Director, Breast Imaging and Intervention
 Gwinnett Breast Center, Gwinnett Medical Center
 Lawrenceville, GA 30046




MammoStar™ at Deployment - Fig. 1



MammoStar™ at 8 Days - Fig. 2



MammoStar™ Wire Localization - Fig. 3



Wire Loc Confirmation - Fig. 4