



Mammotome® Markers

Ordering Chart



Barrel



Butterfly



Open Coil



Bowtie



Triple Twist



U-shape



Barbell



Tribell

Catalog Number	Used With	Gauge	Carrier Plug	Marker Properties	Shape	Ultrasound Visibility	(Modality) & Product Description	
4010-02-18-T3	Mammotome Elite® Biopsy Device & Core Needle	18G	(Mini) PEG Hydrogel	Titanium	Open Coil	Up to 12 months	HydroMARK: rigid 18G needle with radiopaque clip within Hydrogel	
4010-02-15-T1		15G	PEG Hydrogel	Titanium	Barrel	Up to 12 months	HydroMARK: rigid 15G needle with radiopaque clip within Hydrogel	
4010-02-15-T3			PEG Hydrogel	Titanium	Open Coil	Up to 12 months	HydroMARK: rigid 15G needle with radiopaque clip within Hydrogel	
4010-02-15-T4			PEG Hydrogel	Titanium	Butterfly	Up to 12 months	HydroMARK: rigid 15G needle with radiopaque clip within Hydrogel	
4010-02-15-S1			PEG Hydrogel	Stainless Steel	Barrel	Up to 12 months	HydroMARK: rigid 15G needle with radiopaque clip within Hydrogel	
4010-02-15-S3			PEG Hydrogel	Stainless Steel	Open Coil	Up to 12 months	HydroMARK rigid 15G needle with radiopaque clip within Hydrogel	
MAM3014		14G	Bovine Collagen	Titanium	Bowtie	4-8 weeks	CorMARK: rigid 14G needle with radiopaque clip within Bovine Collagen	
MAM3015			Bovine Collagen	Titanium	U-shape	4-8 weeks	CorMARK: rigid 14G needle with radiopaque clip within Bovine Collagen	
STAR1401			Beta Glucan	Ceramic/Carbon	Barbell	Long-term	MammoSTAR: rigid 14G needle with radiopaque clip within Beta Glucan	
STAR1402			Beta Glucan	Ceramic/Carbon	Tribell	Long-term	MammoSTAR: rigid 14G needle with radiopaque clip within Beta Glucan	
STAR1403	(Petite) Beta Glucan		Ceramic/Carbon	Barbell	Long-term	MammoSTAR: rigid 14G needle with radiopaque clip within Beta Glucan		
STAR1001	10G	Beta Glucan	Ceramic/Carbon	Barbell	Long-term	MammoSTAR: rigid blunt tip applicator with radiopaque clip within Beta Glucan		
MMK0801	Mammotome Revolve™ Biopsy Device	8G	Bovine Collagen	Titanium	Bowtie	4-8 weeks	(U/S, ST) MammoMARK: 8G flexible, side-deploy applicator with radiopaque clip within Bovine Collagen	
MMK0802			Bovine Collagen	Titanium	U-shape	4-8 weeks	(U/S, ST) MammoMARK: 8G flexible, side-deploy applicator with radiopaque clip within Bovine Collagen	
MMK0803			Bovine Collagen	Titanium	Triple Twist	4-8 weeks	(U/S, ST) MammoMARK: 8G flexible, side-deploy applicator with radiopaque clip within Bovine Collagen	
STAR0831			Beta Glucan	Ceramic/Carbon	Barbell	Long-term	(U/S, ST) MammoSTAR: 8G flexible, side-deploy applicator with radiopaque clip within Beta Glucan	
STAR0832			Beta Glucan	Ceramic/Carbon	Tribell	Long-term	(U/S, ST) MammoSTAR: 8G flexible, side-deploy applicator with radiopaque clip within Beta Glucan	
STAR0833			(Petite) Beta Glucan	Ceramic/Carbon	Barbell	Long-term	(U/S, ST) MammoSTAR: 8G flexible, side-deploy applicator with radiopaque clip within Beta Glucan	
4010-05-08-T1			PEG Hydrogel	Titanium	Barrel	Up to 12 months	(U/S, ST) HydroMARK: 8G flexible, side-deploy applicator with radiopaque clip within Hydrogel	
4010-05-08-T3			PEG Hydrogel	Titanium	Open Coil	Up to 12 months	(U/S, ST) HydroMARK: 8G flexible, side-deploy applicator with radiopaque clip within Hydrogel	
4010-05-08-T4			PEG Hydrogel	Titanium	Butterfly	Up to 12 months	(U/S, ST) HydroMARK: 8G flexible, side-deploy applicator with radiopaque clip within Hydrogel	
4010-05-10-T1			10G	PEG Hydrogel	Titanium	Barrel	Up to 12 months	(U/S, ST) HydroMARK: 10G flexible, side-deploy applicator with radiopaque clip within Hydrogel
4010-05-10-T3	PEG Hydrogel	Titanium		Open Coil	Up to 12 months	(U/S, ST) HydroMARK: 10G flexible, side-deploy applicator with radiopaque clip within Hydrogel		
4010-05-10-T4	PEG Hydrogel	Titanium		Butterfly	Up to 12 months	(U/S, ST) HydroMARK: 10G flexible, side-deploy applicator with radiopaque clip within Hydrogel		
MMK1001	Bovine Collagen	Titanium		Bowtie	4-8 weeks	(U/S, ST) MammoMARK: 10G flexible, side-deploy applicator with radiopaque clip within Bovine Collagen		
MMK1002	Bovine Collagen	Titanium		U-shape	4-8 weeks	(U/S, ST) MammoMARK: 10G flexible, side-deploy applicator with radiopaque clip within Bovine Collagen		
MMK1003	Bovine Collagen	Titanium		Triple Twist	4-8 weeks	(U/S, ST) MammoMARK: 10G flexible, side-deploy applicator with radiopaque clip within Bovine Collagen		
STAR1031	Beta Glucan	Ceramic/Carbon		Barbell	Long-term	(U/S, ST) MammoSTAR: 10G flexible, side-deploy applicator with radiopaque clip within Beta Glucan		
STAR1032	Beta Glucan	Ceramic/Carbon		Tribell	Long-term	(U/S, ST) MammoSTAR: 10G flexible, side-deploy applicator with radiopaque clip within Beta Glucan		
STAR1033	(Petite) Beta Glucan	Ceramic/Carbon		Barbell	Long-term	(U/S, ST) MammoSTAR: 10G flexible, side-deploy applicator with radiopaque clip within Beta Glucan		
4010-01-08-T1	Mammotome® Biopsy Device	8G		PEG Hydrogel	Titanium	Barrel	Up to 12 months	(U/S, ST, MRI) HydroMARK: 8G flexible, side-deploy applicator with radiopaque clip within Hydrogel
4010-01-08-T3			PEG Hydrogel	Titanium	Open Coil	Up to 12 months	(U/S, ST, MRI) HydroMARK: 8G flexible, side-deploy applicator with radiopaque clip within Hydrogel	
4010-01-08-T4			PEG Hydrogel	Titanium	Butterfly	Up to 12 months	(U/S, ST, MRI) HydroMARK: 8G flexible, side-deploy applicator with radiopaque clip within Hydrogel	
4010-01-08-S1			PEG Hydrogel	Stainless Steel	Barrel	Up to 12 months	(U/S, ST, MRI) HydroMARK: 8G flexible, side-deploy applicator with radiopaque clip within Hydrogel	
MAM3008			Bovine Collagen	Titanium	U-shape	4-8 weeks	(U/S, ST) MammoMARK: 8G flexible, side-deploy applicator with radiopaque clip within Bovine Collagen	
MRM4008			Bovine Collagen	Titanium	Bowtie	4-8 weeks	(U/S, ST, MRI) MammoMARK: 8G flexible, side-deploy applicator with radiopaque clip within Bovine Collagen	
STAR0801			Beta Glucan	Ceramic/Carbon	Barbell	Long-term	(U/S, ST) MammoSTAR: 8G flexible, side-deploy applicator with radiopaque clip within Beta Glucan	
STAR0802			Beta Glucan	Ceramic/Carbon	Tribell	Long-term	(U/S, ST) MammoSTAR: 8G flexible, side-deploy applicator with radiopaque clip within Beta Glucan	
4010-01-11-T1			11G	PEG Hydrogel	Titanium	Barrel	Up to 12 months	(U/S, ST, MRI) HydroMARK: 11G flexible, side-deploy applicator with radiopaque clip within Hydrogel
4010-01-11-T3				PEG Hydrogel	Titanium	Open Coil	Up to 12 months	(U/S, ST, MRI) HydroMARK: 11G flexible, side-deploy applicator with radiopaque clip within Hydrogel
4010-01-11-T4	PEG Hydrogel	Titanium		Butterfly	Up to 12 months	(U/S, ST, MRI) HydroMARK: 11G flexible, side-deploy applicator with radiopaque clip within Hydrogel		
MAM3001	Bovine Collagen	Titanium		Bowtie	4-8 weeks	(U/S, ST, MRI) MammoMARK: 11G flexible, side-deploy applicator with radiopaque clip within Bovine Collagen		
MAM3002	Bovine Collagen	Titanium		U-shape	4-8 weeks	(U/S, ST, MRI) MammoMARK: 11G flexible, side-deploy applicator with radiopaque clip within Bovine Collagen		
MRM4002	Bovine Collagen	Titanium		Bowtie	4-8 weeks	(U/S, ST, MRI) MammoMARK: 11G flexible, side-deploy applicator with radiopaque clip within Bovine Collagen		
STAR1101	Beta Glucan	Ceramic/Carbon		Barbell	Long-term	(U/S, ST) MammoSTAR: 11G flexible, side-deploy applicator with radiopaque clip within Beta Glucan		
STAR1102	Beta Glucan	Ceramic/Carbon		Tribell	Long-term	(U/S, ST) MammoSTAR: 11G flexible, side-deploy applicator with radiopaque clip within Beta Glucan		



Mammotome® Markers


Ordering Chart



Barrel



Butterfly



Open Coil



Bowtie



Triple Twist



U-shape



Barbell



Tribell

Catalog Number	Used With	Gauge	Carrier Plug	Marker Properties	Shape	Ultrasound Visibility	(Modality) & Product Description	
4010-02-18-T3	Core Needle	18G	(Mini) PEG Hydrogel	Titanium	Open Coil	Up to 12 months	HydroMARK: rigid 18G needle with radiopaque clip within Hydrogel	
4010-02-15-T1		15G	PEG Hydrogel	Titanium	Barrel	Up to 12 months	HydroMARK: rigid 15G needle with radiopaque clip within Hydrogel	
4010-02-15-T3			PEG Hydrogel	Titanium	Open Coil	Up to 12 months	HydroMARK: rigid 15G needle with radiopaque clip within Hydrogel	
4010-02-15-T4			PEG Hydrogel	Titanium	Butterfly	Up to 12 months	HydroMARK: rigid 15G needle with radiopaque clip within Hydrogel	
4010-02-15-S1			PEG Hydrogel	Stainless Steel	Barrel	Up to 12 months	HydroMARK: rigid 15G needle with radiopaque clip within Hydrogel	
4010-02-15-S3			PEG Hydrogel	Stainless Steel	Open Coil	Up to 12 months	HydroMARK: rigid 15G needle with radiopaque clip within Hydrogel	
MAM3014			14G	Bovine Collagen	Titanium	Bowtie	4-8 weeks	CorMARK: rigid 14G needle with radiopaque clip within Bovine Collagen
MAM3015		Bovine Collagen		Titanium	U-shape	4-8 weeks	CorMARK: rigid 14G needle with radiopaque clip within Bovine Collagen	
STAR1401		Beta Glucan		Ceramic/Carbon	Barbell	Long-term	MammoSTAR: rigid 14G needle with radiopaque clip within Beta Glucan	
STAR1402		Beta Glucan		Ceramic/Carbon	Tribell	Long-term	MammoSTAR: rigid 14G needle with radiopaque clip within Beta Glucan	
STAR1403		(Petite) Beta Glucan		Ceramic/Carbon	Barbell	Long-term	MammoSTAR: rigid 14G needle with radiopaque clip within Beta Glucan	
4010-04-09-T1	Hologic Eviva®	9G & 12G	PEG Hydrogel	Titanium	Barrel	Up to 12 months	(ST) HydroMARK: rigid blunt tip applicator (with three adapters) with radiopaque clip within Hydrogel	
4010-04-09-T3			PEG Hydrogel	Titanium	Open Coil	Up to 12 months	(ST) HydroMARK: rigid blunt tip applicator (with three adapters) with radiopaque clip within Hydrogel	
4010-04-09-T4			PEG Hydrogel	Titanium	Butterfly	Up to 12 months	(ST) HydroMARK: rigid blunt tip applicator (with three adapters) with radiopaque clip within Hydrogel	
4010-04-09-S3			PEG Hydrogel	Stainless Steel	Open Coil	Up to 12 months	(ST) HydroMARK: rigid blunt tip applicator (with three adapters) with radiopaque clip within Hydrogel	
MAM3001		Bovine Collagen	Titanium	Bowtie	4-8 weeks	(ST) MammoMARK: 11G flexible, side-deploy applicator with radiopaque clip within Bovine Collagen, 20mm aperture - 13cm needle length		
MAM3002						(ST) MammoMARK: 11G flexible, side-deploy applicator with radiopaque clip within Bovine Collagen, 20mm aperture 13cm needle length		
MAM3008		Bovine Collagen	Titanium	U-shape	4-8 weeks	(ST) MammoMARK: 11G flexible, side-deploy applicator with radiopaque clip within Bovine Collagen, 12mm aperture 13cm needle length		
STAR1121		Beta Glucan	Ceramic/Carbon	Barbell	4-8 weeks	(ST) MammoSTAR: flexible, end-deploy applicator with radiopaque clip within Beta Glucan, 20mm aperture 13cm needle length		
4010-01-11-T1		Hologic ATEC® (detached needle)	9G	PEG Hydrogel	Titanium	Barrel	Up to 12 months	(U/S, ST) HydroMARK: 11G flexible, side-deploy applicator with radiopaque clip within Hydrogel
4010-01-11-T3				PEG Hydrogel	Titanium	Open Coil	Up to 12 months	(U/S, ST) HydroMARK: 11G flexible, side-deploy applicator with radiopaque clip within Hydrogel
4010-01-11-T4	PEG Hydrogel			Titanium	Butterfly	Up to 12 months	(U/S, ST) HydroMARK: 11G flexible, side-deploy applicator with radiopaque clip within Hydrogel	
MAM3001	Bovine Collagen			Titanium	Bowtie	4-8 weeks	(ST) MammoMARK: 11G flexible, side-deploy applicator with radiopaque clip within Bovine Collagen, 20mm aperture	
MAM3002	Bovine Collagen			Titanium	U-shape	4-8 weeks	(ST) MammoMARK: 11G flexible, side-deploy applicator with radiopaque clip within Bovine Collagen, 20mm aperture	
STAR1101	Beta Glucan			Ceramic/Carbon	Barbell	Long-term	(U/S, ST) MammoSTAR: 11G flexible, side-deploy applicator with radiopaque clip within Beta Glucan, 20mm aperture	
STAR1102	Beta Glucan		Ceramic/Carbon	Tribell	Long-term	(U/S, ST) MammoSTAR: 11G flexible, side-deploy applicator with radiopaque clip within Beta Glucan, 20mm aperture		
4010-03-15-T1-SHORT	12G		PEG Hydrogel	Titanium	Barrel	Up to 12 months	(U/S, ST) HydroMARK: rigid, side-deploy applicator with radiopaque clip within Hydrogel	
4010-03-15-T3-SHORT			PEG Hydrogel	Titanium	Open Coil	Up to 12 months	(U/S, ST) HydroMARK: rigid, side-deploy applicator with radiopaque clip within Hydrogel	
4010-05-10-T1	Hologic ATEC® (back of probe)		9G	PEG Hydrogel	Titanium	Barrel	Up to 12 months	(U/S, ST) HydroMARK: 10G flexible, side-deploy applicator with radiopaque clip within Hydrogel
4010-05-10-T3		PEG Hydrogel		Titanium	Open Coil	Up to 12 months	(U/S, ST) HydroMARK: 10G flexible, side-deploy applicator with radiopaque clip within Hydrogel	
4010-05-10-T4		PEG Hydrogel		Titanium	Butterfly	Up to 12 months	(U/S, ST) HydroMARK: 10G flexible, side-deploy applicator with radiopaque clip within Hydrogel	
MMK1001		Bovine Collagen		Titanium	Bowtie	4-8 weeks	(U/S, ST) MammoMARK: 10G flexible, side-deploy applicator with radiopaque clip within Bovine Collagen, 20mm aperture	
MMK1002		Bovine Collagen		Titanium	U-shape	4-8 weeks	(U/S, ST) MammoMARK: 10G flexible, side-deploy applicator with radiopaque clip within Bovine Collagen, 20mm Aperture	
MMK1003		Bovine Collagen		Titanium	Triple Twist	4-8 weeks	(U/S, ST) MammoMARK: 10G flexible, side-deploy applicator with radiopaque clip within Bovine Collagen, 20mm aperture	
STAR1031		Beta Glucan	Ceramic/Carbon	Barbell	4-8 weeks	(U/S, ST) MammoSTAR: 10G flexible, side-deploy applicator with radiopaque clip within Beta Glucan, 20mm aperture		
STAR1032		Beta Glucan	Ceramic/Carbon	Tribell	4-8 weeks	(U/S, ST) MammoSTAR: 10G flexible, side-deploy applicator with radiopaque clip within Beta Glucan, 20mm aperture		
STAR1033		(Petite) Beta Glucan	Ceramic/Carbon	Barbell	4-8 weeks	(U/S, ST) MammoSTAR: 10G flexible, side-deploy applicator with radiopaque clip within Beta Glucan, 12mm aperture		
4010-01-08-T1		Hologic ATEC® (MRI Introducer Sheath)	9G & 12G	PEG Hydrogel	Titanium	Barrel	Up to 12 months	(MRI) HydroMARK: 8G flexible, side-deploy applicator with radiopaque clip within Hydrogel
4010-01-08-T3	PEG Hydrogel			Titanium	Open Coil	Up to 12 months	(MRI) HydroMARK: 8G flexible, side-deploy applicator with radiopaque clip within Hydrogel	
4010-01-08-T4	PEG Hydrogel			Titanium	Butterfly	Up to 12 months	(MRI) HydroMARK: 8G flexible, side-deploy applicator with radiopaque clip within Hydrogel	
4010-01-08-S1	PEG Hydrogel			Stainless Steel	Barrel	Up to 12 months	(MRI) HydroMARK: 8G flexible, side-deploy applicator with radiopaque clip within Hydrogel	
MAM3001	Bovine Collagen			Titanium	Bowtie	4-8 weeks	(MRI) MammoMARK: 11G flexible, side-deploy applicator with radiopaque clip within Bovine Collagen, 20mm aperture	
MAM3002	Bovine Collagen			Titanium	U-shape	4-8 weeks	(MRI) MammoMARK: 11G flexible, side-deploy applicator with radiopaque clip within Bovine Collagen, 20mm aperture	



Mammotome® Markers

Ordering Chart



Barrel



Butterfly



Open Coil



Bowtie



Triple Twist



U-shape



Barbell



Tribell

Catalog Number	Used With	Gauge	Carrier Plug	Marker Properties	Shape	Ultrasound Visibility	(Modality) & Product Description		
4010-02-18-T3	Core Needle	18G	(Mini) PEG Hydrogel	Titanium	Open Coil	Up to 12 months	HydroMARK: rigid 18G needle with radiopaque clip within Hydrogel		
4010-02-15-T1		15G	PEG Hydrogel	Titanium	Barrel	Up to 12 months	HydroMARK: rigid 15G needle with radiopaque clip within Hydrogel		
4010-02-15-T3			PEG Hydrogel	Titanium	Open Coil	Up to 12 months	HydroMARK: rigid 15G needle with radiopaque clip within Hydrogel		
4010-02-15-T4			PEG Hydrogel	Titanium	Butterfly	Up to 12 months	HydroMARK: rigid 15G needle with radiopaque clip within Hydrogel		
4010-02-15-S1			PEG Hydrogel	Stainless Steel	Barrel	Up to 12 months	HydroMARK: rigid 15G needle with radiopaque clip within Hydrogel		
4010-02-15-S3			PEG Hydrogel	Stainless Steel	Open Coil	Up to 12 months	HydroMARK: rigid 15G needle with radiopaque clip within Hydrogel		
MAM3014			14G	Bovine Collagen	Titanium	Bowtie	4-8 weeks	CorMARK: rigid 14G needle with radiopaque clip within Bovine Collagen	
MAM3015		Bovine Collagen		Titanium	U-shape	4-8 weeks	CorMARK: rigid 14G needle with radiopaque clip within Bovine Collagen		
STAR1401		Beta Glucan		Ceramic/Carbon	Barbell	Long-term	MammoSTAR: rigid 14G needle with radiopaque clip within Beta Glucan		
STAR1402		Beta Glucan		Ceramic/Carbon	Tribell	Long-term	MammoSTAR: rigid 14G needle with radiopaque clip within Beta Glucan		
STAR1403	(Petite) Beta Glucan	Ceramic/Carbon		Barbell	Long-term	MammoSTAR: rigid 14G needle with radiopaque clip within Beta Glucan			
MRM4008	BD EnCor Enspire™	7G	Bovine Collagen	Titanium	Bowtie	4-8 weeks	(U/S, ST) MammoMARK: 8G flexible, side-deploy applicator with radiopaque clip within Bovine Collagen		
STAR0831			Beta Glucan	Ceramic/Carbon	Barbell	Long-term	(U/S, ST) MammoSTAR: 8G flexible, side-deploy applicator with radiopaque clip within Beta Glucan		
STAR0832			Beta Glucan	Ceramic/Carbon	Tribell	Long-term	(U/S, ST) MammoSTAR: 8G flexible, side-deploy applicator with radiopaque clip within Beta Glucan		
STAR0833			(Petite) Beta Glucan	Ceramic/Carbon	Barbell	Long-term	(U/S, ST) MammoSTAR: 8G flexible, side-deploy applicator with radiopaque clip within Beta Glucan		
4010-05-08-T1			PEG Hydrogel	Titanium	Barrel	Up to 12 months	(U/S, ST) HydroMARK: 8G flexible, side-deploy applicator with radiopaque clip within Hydrogel		
4010-05-08-T3			PEG Hydrogel	Titanium	Open Coil	Up to 12 months	(U/S, ST) HydroMARK: 8G flexible, side-deploy applicator with radiopaque clip within Hydrogel		
4010-05-08-T4			PEG Hydrogel	Titanium	Butterfly	Up to 12 months	(U/S, ST) HydroMARK: 8G flexible, side-deploy applicator with radiopaque clip within Hydrogel		
4010-05-10-T1		10G	PEG Hydrogel	Titanium	Barrel	Up to 12 months	(U/S, ST) HydroMARK: 10G flexible, side-deploy applicator with radiopaque clip within Hydrogel		
4010-05-10-T3			PEG Hydrogel	Titanium	Open Coil	Up to 12 months	(U/S, ST) HydroMARK: 10G flexible, side-deploy applicator with radiopaque clip within Hydrogel		
4010-05-10-T4			PEG Hydrogel	Titanium	Butterfly	Up to 12 months	(U/S, ST) HydroMARK: 10G flexible, side-deploy applicator with radiopaque clip within Hydrogel		
4010-03-09-T1			PEG Hydrogel	Titanium	Barrel	Up to 12 months	(U/S, ST) HydroMARK: 10G rigid, side-deploy applicator with radiopaque clip within Hydrogel		
4010-03-09-T3			PEG Hydrogel	Titanium	Open Coil	Up to 12 months	(U/S, ST) HydroMARK: 10G rigid, side-deploy applicator with radiopaque clip within Hydrogel		
MRM4002			Bovine Collagen	Titanium	Bowtie	Up to 12 months	(U/S, ST) MammoMARK: 11G flexible, side-deploy applicator with radiopaque clip within Bovine Collagen		
STAR1031			Beta Glucan	Ceramic/Carbon	Barbell	Long-term	(U/S, ST) MammoSTAR: 10G flexible, side-deploy applicator with radiopaque clip within Beta Glucan		
STAR1032			Beta Glucan	Ceramic/Carbon	Tribell	Long-term	(U/S, ST) MammoSTAR: 10G flexible, side-deploy applicator with radiopaque clip within Beta Glucan		
STAR1033			(Petite) Beta Glucan	Ceramic/Carbon	Barbell	Long-term	(U/S, ST) MammoSTAR: 10G flexible, side-deploy applicator with radiopaque clip within Beta Glucan		
4010-01-08-T1			BD EnCor Enspire™ (MRI Introducer)	7G & 10G	PEG Hydrogel	Titanium	Barrel	Up to 12 months	(MRI) HydroMARK: 8G flexible, side-deploy applicator with radiopaque clip within Hydrogel
4010-01-08-T3					PEG Hydrogel	Titanium	Open Coil	Up to 12 months	(MRI) HydroMARK: 8G flexible, side-deploy applicator with radiopaque clip within Hydrogel
4010-01-08-T4					PEG Hydrogel	Titanium	Butterfly	Up to 12 months	(MRI) HydroMARK: 8G flexible, side-deploy applicator with radiopaque clip within Hydrogel
4010-01-08-S1	PEG Hydrogel	Stainless Steel			Barrel	Up to 12 months	(MRI) HydroMARK: 8G flexible, side-deploy applicator with radiopaque clip within Hydrogel		
4010-01-08-T1	BD VACORA® (MRI Coaxial Cannula)	10G	PEG Hydrogel	Titanium	Barrel	Up to 12 months	(MRI) HydroMARK: 8G flexible, side-deploy applicator with radiopaque clip within Hydrogel, 118mm MRI Coaxial Cannula		
4010-01-08-T3			PEG Hydrogel	Titanium	Open Coil	Up to 12 months	(MRI) HydroMARK: 8G flexible, side-deploy applicator with radiopaque clip within Hydrogel, 118mm MRI Coaxial Cannula		
4010-01-08-T4			PEG Hydrogel	Titanium	Butterfly	Up to 12 months	(MRI) HydroMARK: 8G flexible, side-deploy applicator with radiopaque clip within Hydrogel, 118mm MRI Coaxial Cannula		
4010-01-08-S1			PEG Hydrogel	Stainless Steel	Barrel	Up to 12 months	(MRI) HydroMARK: 8G flexible, side-deploy applicator with radiopaque clip within Hydrogel, 118mm MRI Coaxial Cannula		