


Mammotome® Markers




HydroMARK™
Breast Biopsy Site Marker



Three Distinct Shapes

BARREL (T1)



OPEN COIL (T3)




BUTTERFLY (T4)




Available in Titanium (nickel free) and
Stainless Steel options.


Highly Visible Under Ultrasound



Immediately after deployment
• Hyperechoic Artifact




+24 hours after deployment
• Anechoic Artifact




6 months after deployment and
neoadjuvant chemotherapy¹
• Anechoic Artifact

1. Sakamoto, N., Fukuma, E., Tsunoda, Y. et al. Breast Cancer 25, 575–582 (2018). <https://doi.org/10.1007/s12282-018-0854-8>


Visibility Under All Modalities



The HydroMARK™ Marker carrier, hydrogel hydrates to provide up to 12 months of ultrasound visibility.¹



Permanent stereotactic, tomosynthesis, and X-ray visibility.¹




Permanent MRI visibility while uniquely distinguishable on a T2 sequence for up to 12 months.²

2. Shah AD, Mehta AK, Talati N, Brem R, Margolies LR. Clin Imaging. 2019 May-Jun;55:196-212. doi: 10.1016/j.clinimag.2019.05.014. Epub 2019 Jun 14. PMID: 31204189.

Pre-Surgical Localization

Localization under ultrasound provides a faster, more comfortable procedure for the patient.³



HydroMARK™ Marker under ultrasound

The visibility of the HydroMARK™ Marker enables consistent localization under ultrasound.³

3. Blumencranz, P.W., Ellis, D. & Barlowe, K. Ann Surg Oncol 21, 3273–3277 (2014). <https://doi.org/10.1245/s10434-014-3917-x>

Compatibility

Core Needle in Ultrasound Procedures

Direct Puncture - 15G & 18G

Vacuum-Assisted Breast Biopsy Probe*

Mammotome	Mammotome Revolve™ Biopsy Device - 8G & 10G
Hologic®	Eviva®/ Brevera® - 9G & 10G ATEC® - 9G, 12G & MRI
BD	EnCor Enspire™, EnCor Ultra™ - 7G, 10G & MRI

*See the Mammotome Marker compatibility brochure for more information

Mammotome

mammotome.com

1-877-926-2666

Product may not be approved or available in your region. Please check with your local Mammotome representative.

© 2021 Devicor Medical Products, Inc. All rights reserved. Devicor, Mammotome, Mammotome Revolve and HydroMARK are registered trademarks of Devicor Medical Products, Inc., in the USA and optionally in other countries. Hologic, Eviva®, Brevera® and ATEC® are registered trademarks of Hologic, Inc. Bard, EnCor Enspire™ and EnCor Ultra™ are the product names, trade names and/or trademarks of C.R.Bard, Inc., and its affiliates including BD. Other logos, product and/or company names might be trademarks of their respective owners. MDM# 200565 Rev 7/21